

## AUTOCLAVE VERTICAL

### DESCRIPTION

Sturdy Construction, Double walled, Inner chamber i.e. boiler made of 16 SWG stainless steel of 304 quality, outer shell is made of thick mild steel duly painted. The boiler and outer shell is provided with air insulation. Lid is made of thick mild steel duly Stainless steel lined and is tightened by Radial locking system. Lid is mounted on a central screw and is provide with best quality joint less neoprene gasket to avoid any risk and trouble free operation. These autoclaves are fitted with paddle lifting device to open/close the lid. Fitted with glass indicator to indicate the water level inside the boiler. These autoclaves are fitted with pressure gauge, steam release cock, safety valves, spring loaded safety valves etc. These autoclaves can be set at any pressure in between 5 to 22 Provided with drain valve. These autoclaves are hydraulically tested upto 40. Provided with ISI mark heating elements which are easily replaceable and heats the water to desited temperature and pressure. Pressure is indicated by pressure gauge. Supplied complete with stainless steel perforated basket, cord plug to work on 220V +/-90% 50 Hz, single phase.



**MODEL: GEN-001**

SIZE (DIA X DEPTH)	ELECTRIC LOAD	CAPACITY
--------------------	---------------	----------

❖ 250 mm x 450 mm	(1.5 KW)	22ltr
❖ 300 mm x 500 mm	(2 KW)	40ltr
❖ 350 mm x 550 mm	(3 KW)	50ltr
❖ 400 mm x 600 mm	(4 KW)	78ltr
❖ 450 mm x 600 mm	(5 KW)	98ltr
❖ 550 mm x 750 mm	(6 KW)	152ltr



**GenNEXT Lab Technologies Pvt. Ltd**

(An ISO 9001:2015 Certified Company)

Office No. 3, CSC-5, 1<sup>st</sup> Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: [sales@gennextlab.com](mailto:sales@gennextlab.com) & [gennextlab@gmail.com](mailto:gennextlab@gmail.com)

Website: [www.gennextlab.com](http://www.gennextlab.com)