

TOC Analyzer

HTY-MC25 TOC analyzer is a patented product developed independently by Tailin, It is used in testing the total organic carbon in water samples, which has high sensitivity and accuracy.

Model: HTY-MC25



Characteristics

- Membrane / conductivity testing technology, oxidant and UV oxidation technology
- Automatic alarm when exceed the limit value
- Can exclude the interference of halides and hydrocarbons
- Memory the last 12 months record automatically, able to look up the record and print the test result any time
- Operation and maintenance is simple, fast, reliable
- Calibration method is simple, the calibration liquid from our company and the customer are both OK.

Technical Specifications

Power supply	220V /50Hz
Rated output	150W
Measuring limit	0.01mg/L
Measuring precision	±5%
Measuring range	0.01mg/L~25.00mg/L
Analysis time	6 min
Response time	within 30 min
Sample temperature	1~95°C
Dimension	48×42×16 cm
Weight	25 kg