

Auto Sampler

The ASE-18 auto sampler realized auto sampling function when measuring TOC of water sample, it is simple and practical.

This machine is composed by 3 parts as following:

1. Auto rotating platform.
2. Up-down structure.
3. Control system.

Model: ASE-18



Characteristics

1. Easy for operation, quickly and reliable.
2. The operator doesn't need to watch when doing testing.
3. Small, convenience, almost no need to maintain.

Working Principle

1. Using the auto turntable with water sample bottle to change sample.
2. Using up-down structure to realize sampling procedure by sampling probe.
3. Using control system to control motor, to realize auto sampling function.

Technical Specifications

Power supply	220V/50 Hz
Dimensions	31cm x 24cm x 33cm
Sampling Qty	18 piece per time
Sampling time options	15 min, 30 min, 45 min, 60min