

## SEROLOGICAL WATER BATH

### DESCRIPTION

Temperature Range is from ambient to 100 deg c with an accuracy of + 1 deg c. double walled, inner made of Stainless Steel, outer made of heavy duty Mild Steel duly powder coated. Properly insulated. Heating elements are tubular type, made of Copper duly Nickle plated. Covered with perforated S.S. tray. Temperature is controlled with thermostat. Temperature control Knob is graduated in degree centigrade. Provided with On/Off switch, indicator ,cord plug, etc. suitable to work on 220V single phase 50 Hz AC.



**MODEL: GEN-030**

### SIZE (DIA & DEPTH)

300 mm x 250 mm x 175 mm	for 2 racks
300 mm x 300 mm x 175 mm	for 4 racks
455 mm x 300 mm x 175 mm	for 6 racks
455 mm x 455 mm x 175 mm	for 8 racks

### OPTIONAL ACCESSORIES

- ❖ Stirrer with 1.20 HP Motor and speed regulator
- ❖ Digital Temperature Indicator cum controller



## GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015 Certified Company)

Office No. 3, CSC-5, 1<sup>st</sup> Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: [sales@gennextlab.com](mailto:sales@gennextlab.com) & [gennextlab@gmail.com](mailto:gennextlab@gmail.com)

Website: [www.gennextlab.com](http://www.gennextlab.com)