

HOT PLATE

FEATURES:

- ❖ Black backlit 120*55 mm LCD screen
- ❖ Electromagnetic balance sensor
- ❖ External/ Internal calibration
- ❖ 4-level shockproof
- ❖ Counting function
- ❖ Supports PC connection. Wireless Bluetooth printing is optional
- ❖ Top made of cast iron, Body made of mild steel with powder coat. Max. size: 8 inch Dia, Precise temperature control with PID controller, LED display for
- ❖ showing actual and set temperature, Temperature range: Ambient + 5°C to 350°C, Resolution 0.1°C.
- ❖ Temperature accuracy $\pm 0.5^{\circ}\text{C}$
- ❖ Overload protection
- ❖ Input voltage: 220 V AC $\pm 10\%$, 50-60 Hz
- ❖ Rating 1-2 kw



MODEL: GEN-067



GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015 Certified Company)

Office No. 3, CSC-5, 1st Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: sales@gennextlab.com & gennextlab@gmail.com

Website: www.gennextlab.com