

VACUUM OVEN

DESCRIPTION

These are double walled units with outer made of Mild Steel sheet Duly powder coated paint & inner made of heavy gauge Stainless Steel sheet. Gap between the two walls is duly filled with high grade Glass wool. Temperature range from ambient to 150°C is controlled By Digital Temperature Controller with an accuracy of $\pm 1^{\circ}\text{C}$. It is Capable of with standing a high vacuum. Provided with vacuum Gauge, See through transparent window of toughened Glass & one Shelf. To work on 220/230 volts AC supply.



MODEL: GEN-007

INSIDE SIZE (DIA x DEPTH)

- ❖ Round: 200 x 250 mm (8"x10")
- ❖ Round: 250 x 300 mm (10"x12")
- ❖ Round: 300 x 300 mm (12"x12")

OPTIONAL ACCESSORIES

- ❖ Microprocessor PID Control



GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015 Certified Company)

Office No. 3, CSC-5, 1st Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: sales@gennextlab.com & gennextlab@gmail.com

Website: www.gennextlab.com