

HEATING MANTLE (MADE FROM GLSSYARN)

DESCRIPTION

It with stands high temperature up to 350 det C. the body of the mantle is spun one piece from non rusting Aluminium (upto 5.0 liter Capacity) duly Stoving painted. All heating mantels are lagged with special grade mineral wool. To work on 220/230 volts AC.

WITH ENERGY REGULATOR



MODEL: GEN-069

CAPACITY

- ❖ 50 ml or 100 ml
- ❖ 250 ml or 500 ml
- ❖ 1 ltr
- ❖ 2 ltr
- ❖ 3 ltr
- ❖ 5 ltr
- ❖ 10 ltr
- ❖ 20 ltr



GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015 Certified Company)

Office No. 3, CSC-5, 1st Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: sales@gennextlab.com & gennextlab@gmail.com

Website: www.gennextlab.com