

KJELDAHL DISTILLATION UNIT

DESCRIPTION:

It is very useful for estimation of Nitrogen content in the substance by digestion for 300 & 500 ml. Capacity flasks provided with LEAD FUME DUCT and pair of hangers. To work on 220/230 AC. Supply.



MODEL: GEN-071

DELUXE MODEL

- ❖ Unit of 3 test Mantle Type (Without glass parts)
- ❖ Unit of 6 test Mantle Type (Without glass parts)



GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015 Certified Company)

Office No. 3, CSC-5, 1st Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: sales@gennextlab.com & gennextlab@gmail.com

Website: www.gennextlab.com