

MEMBRANE FILTRATION ASSEMBLY

DESCRIPTION:

Complete assembly consists of filter glass Funnel with S.S. Clamp & neoprene lid & neoprene cork. Suitable For Filter Dia 47 mm testing through vacuum filtration. Supplied Without vacuum Flask Cap. 25 Ltr.



MODEL: GEN-064/92

COMPLETE ASSEMBLY CONSISTS OF

- | | |
|---|----------|
| ❖ Filter Funnel S.S.316 with neoprene cover capacity 300ml | one |
| ❖ Membrane Filter C.N. 47mm/0.45um pack of 100 | one pack |
| ❖ Vacuum Flask one liter Borosil | one |
| ❖ Vacuum pump diaphragm type oil free | one |
| ❖ Air Vent Filter Hydrophobic PTFE 47mm/0.2um | one |
| ❖ Syringe type Nylon Membrane Filter 25mm/0.45 um with needle | one |
| ❖ Flexible Tubing 1Mtr. | one |
| ❖ Wooden Box to house all the above | one |



GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015 Certified Company)

Office No. 3, CSC-5, 1st Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: sales@gennextlab.com & gennextlab@gmail.com

Website: www.gennextlab.com