

FLOCCULATOR (JAR TESTING APPARATUS)

DESCRIPTION

Estimate the dosing of alums and such other coagulants for treatment of water and effluent. samples is easily viewed by diffused base illumination with florescent tube placed below plastic plate & black background. Fitted with D.C. motor of 1/20 HP with variable speeds up to 180 RPM controlled by speed regulator. Provided with stainless steel paddles with adjustable spacers to adjust the depth of paddles.

- ❖ For two stirrers
- ❖ For four Stirrers
- ❖ For six Stirrers



MODEL: GEN-050

OPTIONAL ACCESSORIES

- ❖ Digital RPM
- ❖ Automatic 0-60 minute" GERMAN" timer



GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015 Certified Company)

Office No. 3, CSC-5, 1st Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: sales@gennextlab.com & gennextlab@gmail.com

Website: www.gennextlab.com