

## ROTARY VACUUM PUMP

### DESCRIPTION

Vacuum pumps are mounted on a heavy steel base plate, and have a wide range of oil sealed pumps. Fitted with heavy duty electric motor, all pumps are v-belt driven at a speed of 1400 RPM. Pumps are fabricated out of high grade cast iron and are worked on sophisticated machine to give accuracy in microns, thus leading to very high sustained vacuum. Supplied in single stage or double models. All vacuum pumps are fitted with belt guard, vacuum gauge as a standard feature to work on 220/230 volt AC supply.



**MODEL: GEN-063**

### CAPACITY

- ❖ 25 ltr/min
- ❖ 50 ltr/min
- ❖ 75 ltr/min
- ❖ 100 ltr/min
- ❖ 200 ltr/min
- ❖ 300 ltr/min

### MOTOR CAPACITY

- 0.25 H.P.
- 0.25 H.P.
- 0.25 H.P.
- 0.50 H.P.
- 1.00 H.P.
- 1.00 H.P.

### WITH PERFECT/ EQUIVALENT MOTOR DOUBLE STAGE SINGLE STAGE

### OPTIONAL ACCESSORIES

- ❖ Moisture Trap.
- ❖ Vacuum gauge with regulator
- ❖ Non return valve
- ❖ Crompton motor in place of perfect motor



## GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015 Certified Company)

Office No. 3, CSC-5, 1<sup>st</sup> Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: [sales@gennextlab.com](mailto:sales@gennextlab.com) & [gennextlab@gmail.com](mailto:gennextlab@gmail.com)

Website: [www.gennextlab.com](http://www.gennextlab.com)