

WATER BATH RECTANGULAR (DOUBLE WALLED)

DESCRIPTION:

Double walled, inner made of 304 stainless steel, Outer made of Mild Steel duly powder coated. Top made of stainless steel with 75 mm thick holes with concentric rings. Insulation with thick glass wool to prevent any heat loss. Temperature range from ambient to 99 deg c with an accuracy of + 1deg c. thermostatically controlled fitted with Neon Lamp, Indicator main switch and workable on 220V.50 Hz single phase AC. provided with drain and port for thermometer.



MODEL: GEN-029

INNER CHAMBER MADE OF STAINLESS STEEL AND AVAILABLE IN FOLLOWING SIZE.

- ❖ 300mm x 250mm x 100mm (6 hole of 75 mm dia)
- ❖ 300mm x 350mm x 100mm (9 hole of 75 mm dia)
- ❖ 350mm x 400mm x 100 mm (12 hole of 75 mm dia)

OPTIONAL ACCESSORIES

- ❖ Stirrer with Motor and speed regulator
- ❖ Digital Temperature indicator cum controller



GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015 Certified Company)

Office No. 3, CSC-5, 1st Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: sales@gennextlab.com & gennextlab@gmail.com

Website: www.gennextlab.com