

RECTANGULAR HOT PLATE (WITH CAST IRON TOP)

DESCRIPTION

For continuous heating up to 350 deg C. the body is made of specially designed perforated P.C.R.C. Sheet painted with attractive stove enamel with insulated elements. To work on 220/230 volts AC.



MODEL: GEN-068

SIZE (B X L X H)

- ❖ 20 x 30 x 15 cm (10"x 12"x 6")
- ❖ 25 x 40 x 15 cm (10"x 16"x 6")
- ❖ 30 x 45 x 15 cm (12"x 18"x 6")
- ❖ 45 x 60 x 15 cm (18"x 24"x 6")

OPTIONAL ACCESSORIES

- ❖ Thermostate controller.
- ❖ Digital Temperature Controller.



GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015 Certified Company)

Office No. 3, CSC-5, 1st Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: sales@gennextlab.com & gennextlab@gmail.com

Website: www.gennextlab.com