

GLASS DRYER

Features

Benchtop /wall hanging heavy duty electrically driven blower Glass Dryers useful for drying laboratory glassware etc. fast , time saving and safely.

Fully controlled through an energy regulator heating coils blower heat the blown air through the tubes to dry the glassware etc.



MODEL: GEN088

Technical Specification

- ❖ Number of PEGS: 21
- ❖ PEG Size: 16mm x 8 Nos, 10mm x 8 Nos, 7mm x 5 Nos or Customized
- ❖ Temperature Range: Ambient to 70°C
- ❖ Time Control: Yes
- ❖ Heater Selector: Manual
- ❖ Drying Power: 650 W
- ❖ Power: 220V A.C, +/- 10V



GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015, MSME & CE Certified Company)

Office No. 3, CSC-5, 1st Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: sales@gennextlab.com & gennextlab@gmail.com

Website: www.gennextlab.com

