

ANALYTICAL WEIGHING BALANCE

High Performance

Fast, Reliable, Stable Results: Smart Grid weighing pan speeds up your weighing process/Applications. The Smart Grid structure allows air turbulence to flow through the pan, reducing settling and stabilization times for faster results.

Features

- ❖ Colorful Large touch screen
 - ❖ Powerful functions like weighing, parts counting, percent weighing, check weighing, animal weighing, density weighing, etc.
 - ❖ Standard anti-static module
 - ❖ Electronic level adjustment indicator module
- GLP/GMP



MODEL: GEN2204

Technical Specification

- ❖ Display: Touch Screen
- ❖ Weighing Range: 220g
- ❖ Calibration: External
- ❖ Temperature: 10 to 45°C
- ❖ Readability: 0.1mg
- ❖ Response time: 2 to 4 Sec
- ❖ Repeatability: $\pm 0.0002g$
- ❖ Linearity: $\pm 0.0003g$
- ❖ Pan Size: $\Phi 90mm$
- ❖ Interface: RS 232
- ❖ Dimension/N.W: 360x220x350mm/8kg
- ❖ Shipping Size/G. W: 520x380x500mm/10kg



GenNEXT Lab Technologies Pvt. Ltd

(An ISO 9001:2015, MSME & CE Certified Company)

Office No. 3, CSC-5, 1st Floor, Gole Market, Sector-9, Rohini, Delhi – 110085

Email: sales@gennextlab.com & gennextlab@gmail.com

Website: www.gennextlab.com

